REMARKS

Please reconsider the application in view of the above amendments and the following remarks. Applicant thanks the Examiner for carefully considering this application and for indicating that claims 11-16, 18-25, and 31 contain allowable subject matter.

Disposition of Claims

Claims 1-34 are pending in this application. Claims 1 and 32 are independent. The remaining claims depend, directly or indirectly, from claims 1 and 32.

Rejection(s) under 35 U.S.C. § 102

Claims 1-3 and 32 are rejected under 35 U.S.C. § 102(a) as being anticipated by US Patent No. 6,885,801 ("Shankar"). This rejection is respectfully traversed.

Applicant notes that in the Advisory Action mailed June 22, 2009, the Examiner acknowledges that the conception date for the present invention occurred before the priority date of the Shankar reference. Applicant thanks the Examiner for considering the revised 37 C.F.R. § 1.31 Declaration ("revised Declaration") previously filed and acknowledging the present invention's earlier conception date.

However, with respect to diligence, the Examiner continues to assert that the exact date of reduction to practice and the evidence showing diligence is lacking. In response, Applicant respectfully points out that the exact date of reduction to practice is the constructive reduction to practice performed by filing the foreign application to which the present application claims priority, on July 18, 2002.

Further, the Examiner asserts that the evidence documents submitted are not translated. In response, Applicant hereby submits translations of relevant portions of five of the ten original

evidence documents submitted on June 2, 2009. A spreadsheet provides the name, date, and title of each document of which portions have been translated. Translations were made by co-inventor Aymeric Perchant. The translated portions correspond directly to material discussing image processing, which is the concept of the present invention.

Further, the Examiner admits in the Advisory Action dated June 22, 2009, that the last page of the laboratory notebook page indicate that the algorithm was still being “tweaked.” Applicants assert that this supports Applicants’ position that diligence continued after August 2001 (when the laboratory notebook pages end) and was carried through until the filing of the foreign patent application on July 18, 2002. As evidence of furtherance of the concept of the present invention, the translated diligence documents show diligence from November 2001 (the earliest date of the documents provided) until March 2002 (the latest date of the documents provided).

The Examiner asserts that inventor Frederic Berier’s contribution is not shown in the revised Declaration. Applicant respectfully asserts that there is no requirement in the rules under MPEP § 2138.06 that each inventor’s contribution be shown when proving diligence. Further, each inventor’s contribution does not have to be equal. *See* MPEP § 2138.06.

In view of the above, the Applicant asserts that all evidence required by the Examiner, including translated pages of relevant laboratory notebook pages clearly showing a date of conception, and diligence evidence, is submitted by this reply. Accordingly, favorable consideration of the evidence documents detailing diligence from November 2001 to August 2002 is respectfully requested.

Rejection(s) under 35 U.S.C. § 103

Claims 4-5 stand rejected under 35 U.S.C. § 103(a) as being unpatentable over Shankar in view of NPL document “Hybrid Image Segmentation Using Watersheds and Fast Region Merging” (hereinafter “Harris”). Claims 7-10 stand rejected under 35 U.S.C. § 103(a) as being unpatentable over Shankar in view of Harris, and further in view of US Patent No. 4,926,257 (“Miyazaki”). Claims 26-30 and 34 are rejected under 35 U.S.C. § 103(a) as being unpatentable over Shankar in view of US Patent No. 5,764,809 (“Nomami”). For the reasons set forth below, these rejections are respectfully traversed.

Each of the aforementioned §103 rejections are based on Shankar as the primary reference. As described above, Shankar is not valid prior art to this application, as evidenced by revised Declaration filed on June 2, 2009, and the translations of relevant portions of the diligence documents. In view of the above, Shankar may not be used to support the aforementioned 35 U.S.C. § 103 rejections. As the various secondary prior art references admittedly fail to teach all of the elements of the claims, as is evidenced by the fact that they are cited solely to teach specific features recited in dependent claims, it is clear that the various secondary prior art references, whether considered alone or in combination, cannot render the claims of the present application obvious.


In view of the above, withdrawal of the aforementioned rejections with respect to the dependent claims is respectfully requested.

Conclusion

Applicant believes this reply is fully responsive to all outstanding issues and places this application in condition for allowance. If this belief is incorrect, or other issues arise, the Examiner is encouraged to contact the undersigned or his associates at the telephone number listed below. Please apply any charges not covered, or any credits, to Deposit Account 50-0591 (Reference Number [17542/005001]).

Dated: September 8, 2009

Respectfully submitted,

By 
Jonathan P. Osha
Registration No.: 33,986
OSHA · LIANG LLP
909 Fannin Street, Suite 3500
Houston, Texas 77010
(713) 228-8600
(713) 228-8778 (Fax)
Attorney for Applicant

Attachments (Translations of diligence documents)

Filename	Month	Date	Title	Authors	Description/link with image processing	Translated pages
MEMO_acquisition NOV2001.ps	11	2001-11-27	Image Acquisition in November 2001	Georges Le Goualher	Image Acquisition protocol with the tomoscope prototype to acquire proper raw image for image processing testing.	Objective, the protocol
MEMO_devlog.ps	2	2002-02-01	Software Development	Georges Le Goualher	Software architecture, including image processing.	Intro, section 6, drawing
ManuelUtilisation V1.2	2	2002-02-04	Cellvizio User Guide V1.2	Georges Le Goualher, Aymeric Perchant	Software guide as a web page. Screenshots show the fiberdetect functionality, the background subtraction.	Application.html: Sec. 3
ReferenceImageCell-2002-02-07.ps	2	2002-02-07	Reference Guide of ImageCell module	Aymeric Perchant	image reconstruction. Results from the Nov. 2001 database and algorithm block-diagram are shown. First release.	Page 5 to 23
FMODIF interface_cellViZio.ps	3	2002-03-18	Functional Requirement of Software Interface for Software prototype	Georges Le Goualher, Aymeric Perchant	Requirements of the first software prototype controlling the tomoscope, former name of Cellvizio prototype device, and including image processing capabilities.	Introduction, page 2, par on Image Processing

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Modification of the graphical interface of the acquisition module for the Cell-vizio

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1. Introduction

This document in the specifications of the modification of the current graphical interface of the Cell-vizio prototype, for the conformance of the prototype with current regulations (APAVE).

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Image Processing:

The image processing window is kept: the user can chose the number of image to use in the buffer, he can launch the background acquisition using the background button, subtract the background using the subtraction button, and launch image processing using ImageCell. However, the evolution of the FiberDetect and TempMean buttons will depend on the technical feasibility of the processing during acquisition.

We should work on a system that could, during acquisition, assess if the fiber detection previously estimated is correct, or if it has to be re-estimated.

For the temporal averaging, the best would be to compute it regularly, and subtract it to the acquired images ; we still have to know if it can be done in real time... If yes, the TempMean button will become a "Toggle button" : the averaging will be subtracted if and only if the button is enabled.

Therefore, the interface will contain more than four processing buttons: Background, Subtraction, TempMean and ImageCell.

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Software for Cell-ViZio (aka tomoscope in project)

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Introduction

The Software associated with Cell-viZio is launched from the icon mktTest.

The software presents a main window (mktApp) with a browser visualization (modBrowser), and can launch the acquisition module (modACQ) and a visualization module with processing (modPROCESSING)

Standard use scenario


1. Open software using the mktTest icon
2. Using button  create a new "patient". Enter his identification. This action creates a new "patient" and initialize the directories \$Name/Exam1/Study1 in the folders.
3. Select the "Study1" in the third column (note that you can rename this icon with a right click on the mouse, or create a new one with the icon ).
4. When a study is selected, the icon  becomes active
5. Click on this icon to launch the acquisition module
6. A dialog window opens : enter the image prefix that you will use for storage, or click OK to use the default prefix.

7. The acquisition module is displayed full screen. Click on the **Start/Stop** button to launch acquisition
8. When the image is relevant, click on **Pause/Continue** (or the left button of the footswitch) to freeze the image and/or continue
9. Click on the **Save** button or on the right button of the footswitch to save the image with its image buffer. Is the **Save** button is used, a dialog box opens: enter the suffix for the image you want to save (e.g. surf) or **OK** for the default suffix (a number). If the footswitch is used, a default suffix (number) is used.
10. After image saving, the acquisition resumes automatically
11. At the end of the session, click on the **Start/Stop** button to stop acquisition
12. Quit the acquisition module with the close button from the main menu to come back to the browser.
13. Double click on the image icon to visualize it in the modProcessing module.

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3. Section "Processing" (for MKT expert only)

"Background": permits the acquisition of the ghost reflexion for subtraction.

"Subtraction": Permits the subtraction of the background (ghost reflection).


"Temp. Mean": Permits the computation of the temporal averaging of the buffered images.

"Gamma LUT": Permits a gamma transform to improve image contrast

"FiberDetect": launch the process of fiber detection and calibration. This button uses the acquired background if the subtraction button is enabled. If the subtraction mode is disabled, the image background is automatically estimated. At the end of the processing the fiber image with the normalized injection rate is displayed. A white fiber indicates a good injection rate, a dark fiber a weak injection.

"ImageCell": this button is enabled if a "fiberdetect" was activated before. The detected fibers are used for image reconstruction of the acquired data, correcting optical vignetting and injection defects. Working like the fiberdetect button, the background is used only if the subtraction mode is enabled. Instead, the background is estimated automatically.

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!cf. display.exe

1. Introduction

The aim of this document is to define the software developments tasks on the main tomoscope application **mktApp**. The italic text points the first priorities. The references M1, M2, etc... are links to the meeting report on images from E. Borotto.

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3. Modification for the acquisition module: modAcq

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5. Modification for the application: mkt App

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Real time subtraction of the background due to the parasite reflections. The acquisition module permits the storage of the background linked to the ghost reflection, and the computation of the temporal averaging. This background is subtracted from the current image (point to point subtraction between the raw image and the temporal averaging of the background). This permits the enhancement of the current image contrast. The user should be trained for the background acquisition. For the moment, the subtraction is only valid when the acquisition parameters of the background and the current image are identical (problem of the electronic offset, unknown gain control curve). **Should we talk about the footer offset ? : influence on the visualization module ?**
Correction of the offset between lines: it is the default mode of operations (electronic parameters group). **Evolution of the new electronics and the new box ?**

- **Temporal averaging:** today, only this type of averaging can be applied. Is it meaningful to let the user apply other types of averaging ?
- Automated gain adaptation:** it can still be applied due to the fluctuations of the signal linked to the electronic offset. One or two frames could be lost during acquisition.
- **Fiber detection, calibration, reconstruction...** : scenario, interface and implementation ? Choice: with or without fiber detection ?

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Fiber detection, calibration, reconstruction: scenario, interface and implementation ?

Correction of geometrical distortions: can we imagine a processing module in which we can input the necessary informations (cf. prototype software) and a distortion grid pattern image. A deformation field is computed (X and Y). This processed field and the necessary informations are used to remove the distortions on the images. We can then apply the specific measurement tools (interactive ones) to the processed image.

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Subject: Reference Manual of the libMKTPProcessing part of ImageCell. This manual contains informations on image processing, C++ software code and usage.

Introduction

This document gathers different pieces of informations for the main software module entitled Image-Cell. This modules contains different parts of the image processing software:

- fiber detection
- estimation of the optical flow coming from the fibers
- compensation of the device artifacts (bias, background,...)
- device calibration
- image reconstruction.
- Those parts are addressed in details in each chapter to explain algorithms and parameters (chapter 1), their implementation (2) and their usage (3).

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Image Processing in Image-Cell

Introduction

This part describes the algorithms beginning with their generic descriptions, until the individual elements.

1.1 General Schematic

Figure 1.1 presents the whole picture of the processing/ We can distinguish 4 blocks, with 2 identical blocks but their parameters. There exists 3 groups of processing:

Fiber Detection: this processing permits the detection and the segmentation of each fiber on an image, and the analysis of the spatial organization of the fiber bundle;

Flow estimation: a raw image, size 640x640x20 (width, height, number of temporal frames), that contains 8 mega-pixels representing the information viewed by 10,000 or 30,000 fibers. This bloc can isolate data seen by each fiber.

Calibration and reconstruction: The data seen by each fiber is then processed to reconstruct an image without artifacts from the device.

Each of the 3 blocks are detailed hereafter. The prototyping of the algorithms described here are described in the internship report of Sandra Marti ; we encourage the reader to read it for further information on the history of the algorithms, especially the fiber detection one.

The flow estimation left block permits the measurement of the flow on an image representing the spatially constant properties of a given object (a homogeneous diffusive medium, or nearly homogeneous object in random motion, or a mirror). The output of this block is an image of the fiber injection rate in each fiber, taking into account the whole chain of processing.

The flow estimation right block permits the measurement of an object observed through the device. The last block calibration is the dividing of the object image by the fiber rate image. This processing allows for the compensation of the heterogeneous injection in some fibers. Then the image is reconstructed.

1.2 Fiber detection block.

The fiber detection is organized around 4 sequential processing:

- preprocessing,


Fig 1.1 Schematic of the processing. Each bold box are inputs. Dashed box are data or optional processing (bias, background subtraction)

- watershed (LPE) : the first region segmentation,
- split of the too big fibers (split)
- merge of the too small fiber (merge).

The last two steps can be iterated ; nevertheless a performance and convergence analysis should be done. We won't address this matter yet.

1.2.1 Pre-processing

The pre-processing are the following ones.

Anisotropic diffusion: the aim is to smooth the image into flat areas, that is inter-fiber areas. The prototype worked and did not used the inimage library but a library Sandra used issued by ENST. The parameters of the diffusion is postponed when the next release of libinimage contain the proper documentation ! Parameter `[_doDiffusion_]`. **Interpolation x2 nearest neighbor**: we want to simulate structuring elements for mathematical morphology with a radius lower than 1. The image size is doubled so that an opening does not modify the isolated maxima, but only the 8-connected and not 4-connected (neighbor on the diagonal). The usefulness lies in the selection of maxima that are not eliminated by the opening. Parameter `[_doZoom_]`.

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Numerical opening: we want to erase spurious maxima located on the fibers. This is a standard pre-processing of watersheds, possibly modified by the behavior of the size doubling previously described.

Image inversion: the watershed works with minima. Therefore, we invert the image to turn maxima into minima.

1.2.2 First watershed

The watershed of libInImage needs the definition of proper landmarks in the image. We obtain landmarks as a minimum of the pre-processed image using h-domes. For minima, h-domes are defined for an integer h as a function f : $Z \rightarrow Z$ for the image:

$$g = E_f(f + h) \quad f,$$

with $E_a(b)$ being the geodesic reconstruction by erosion of b in a . g is the set of h-domes, that we threshold to one to retrieve all of them. We obtain the minima of depth at least equal to h . This h parameter is given by `[hDomeWatershed_]`.

The result of the watershed is a connected component image containing the detected fibers. To avoid border problems, we glue component touching the borders to the borders. Those components are with the borders are set to the connected component with index 0, the reserved index for the background. Last, the connected components are re-organized to fill the holes, and are sorted with their sizes.

1.2.3 Split

The pre-processing, and the numerical opening are used to limit the main drawback of the watershed: oversegmentation. The effect is that several segments will gather several objects (the fibers). The first defect that we compensate for is the inevitable under-segmentation due to the processing, and even added by the pre-processing. The processing achieved here are:

Selection of under-segmented suspected segment: we found the fibers the size of which is greater than a given threshold depending on the average normalized size (parameter `[_oversizedFactor_]`, and depending on the number of neighbors `[_overSizedNeighLimit_]`).

Re-segmentation: the selected segments are isolated and re-segmented either on the source image without pre-processing, or on a distance map inside those regions (parameter `[_overSizedResegmentMethod_]`). In the first case, we suppose that the main defects come from the pre-processing, in the second case we suppose that they were here before, and that we want to split the segment based on the narrowing of the shape (shape of a peanut). The watershed is done on the same parameters as the first one.

1.2.4 Merge

Then we want to compensate for the main defect (for a watershed): oversegmentation. We want to merge segments. To do that, we do:

- preselection of merging candidates
- among them, find the ones that will be segmented (mandatory)
- among the other ones : eliminate the impossible fusions, and list the possible fusions.

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- éliminer la merging with poor results
- if there is no more merging after filtering, pull out the fiber

For all other fibers, take the best merging that results in the best compactness for the merged segments.

The important points are detailed hereafter:

Remarks on compactness

The compactness that is used is the most common. For an object with perimeter P , and surface S , the criteria is P^2 / S . The perimeter is computed as the 4-connected border, et is defined for all pixel objects having an 8-connected neighborhood with the background (defined here as the rest of the image). The surface is the object itself.

However, this simple criteria is not invariant to scale. We can show that for simple shapes, the criteria is increasing with size, instead of being constant. But the average increase in only in $1/\sqrt{n}$, with n being the diameter of the object (here, for a sphere in 4-connectivity). This criteria is still reasonable, but more subtle to deal with when comparing object of various sizes.

One have to note that for the first trials on a distorted image we should either compensate for them, or take into account the spatially variant anisotropy : we would over-estimate both the border and the surface on the borders of the image.

Selection of the too small detected fibers

The selection has 2 steps. The first one is the selection of the potential segments for merging using a threshold on size. We use the parameter `[_underSizedFactor_]` that is a multiplicative factor of the mean size, and that is generally lower than 1. The second step is the selection among those segments of the ones that will be merge, et the ones that won't.

Mandatorily merged segments

The selection of the segments (or fibers) that will mandatorily be merged is done in the space (number of neighbors) x (normalized size). The size is normalized against the average (we divide all the sizes by the average size). The number of neighbors is computed using 8-connectivity.

Figures 1.2 and 1.3 show the joint histogram of the normalized size and the number of neighbors for each fiber. We remark that their exists a main axis centered on the the line linking the point of interest (6, 1), (8, 2), and (5, 5). Those points can be explained geometrically as indicated on figure 1.4.

This axis is a natural axis for the hexagonal structure of the fibers in this space. If we orient this space, the little fibers will be located towards a lower number of neighbors, and towards smaller fibers. We fix a linear border to select the too small detected fibers. After a few tests on real images of the November 2001 database, and depending on the preliminary results of [?], we choose a line defined by the points (6, $1/3$) (3, 1). The parameters of this line are given by the equation $y = mx + p$ with m the slope (`[_underSizedMandatoryMergeSlope_]` and p the intercept (`[_underSizedMandatoryMergeOffset_]`). All points below this line will be fused, and won't go through the filters described hereafter.

Filtering the potential merging

We define filtering here a process to remove possible merging. There are 3 possible filters.

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Fig 1.2: 3D view of the 2D joint histogram of the normalized size et the number of neighbors for each fiber. The geodesics are projected below: the lines are centered around 6 neighbors, and the average size.

Filtering on the merged size: (parameter `[_underSizedFilterSize_]`) we remove the merging the size of which is above a given threshold fixed by `[_underSizedFilterSizeMax_]`, and that is a multiplicative factor applied to the average size.

Filtering on the total compacity: (parameter `[_underSizedFilterCompact_]`) we remove the merging the compacity of which is above a given threshold fixed by `[_underSizedFilterCompactMax_]`, that is a multiplicative factor of the standard deviation added to the average size.

Filtering on the change of compacity: (parameter `[_underSizedFilterMoreCompact_]`) we remove the merging that increases the compacity of the merged object.

For the default parameters, we use only the first filtering that seems enough in most cases. The other filters seem to be too restrictive et should be studied more.

1.3 Flow estimation block

The flow estimation block has subparts:

Flow estimation as seen by each fiber: this is what we want to retrieve from an image using the fiber detection;

Estimation of the image background: this background is the ghost reflections or (including!) the electronic offset of the detector;

Background subtraction: we remove the background of the image

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
Fig 1.3: the geodesics of the 2D joint histogram of the normalized size and the number of neighbors for each fibers. The 3d shape of this histogram is presented on figure 1.2


Fig 1.4 - Variability of the number of neighbors depending on the size of the detected fiber.

Bias correction: images are often biased, that is there exists a slowly varying background, and not a constant background.

The subtraction block is optional, and also the bias correction bloc,

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1.3.1 Flow estimation for each fiber

The report [?] describes a method to compensate for saturation issues on images. Those problems will be solved by the new electronics arriving mid may, 2002. Therefore we have decided to stop investigating this path for now. The flow is estimated using only an averaging.

Several enhancement are possible:

- consider the Poisson noise in the max estimator using an averaging,
 - consider the fiber shape disparity of the fiber bundle: all fibers are different, et the flow estimation should be adapted to each fiber,
 - check possible issues due to the inter-core spacing, depending on the focal spot on the fiber bundle,
 - more generic: find the complete model for the injection and the backscattering of the laser spot described in [?].
- Those enhancements could be studied, but have low priority for now.

1.3.2 Background estimation or usage, then subtraction

The background can have several sources: they are described in the memos [?, ?, ?]. To sum up, we call background here all ghost reflexions on optics, including the fiber bundle distal end, and the electronic offset of the digital converter.

If the offset is the biggest part on the image, we cannot measure the background simply by subtracting the image because the offset depends on the content, and differs. In this case, we use a percentile of the histogram ($1/100$ or $1/1000$) for its estimation.

Otherwise, we need to use the acquired background by removing the observed object, and then subtract it. In all cases, we need to threshold the subtraction to avoid any artifacts due to negative outliers. The use of a footer (cf [?]) is not maintained because the diffusion background is not a main component of the signal.

1.3.3 Bias correction

The bias should *a priori* be compensated on the injection rate acquisition (left hand box) and on the object acquisition (right hand box). In the first case, this comes from the fact that the calibration is done on a plan mirror, and that the field curvature (??) will reduce the injection quality on the back scattering on the borders (not focalized). This should not be the case.

On the object, or in any homogeneous diffusive medium, the injection stay lower on the borders, and this result in a bias very similar to the first one.

In all cases, the bias as a circular symmetry.

The estimation of the bias is done by dividing the image in $N \times N$ blocks with the same size, then estimating the bias on each block. We need to consider the object nature. On a homogeneous object, the bias can be acquired using the mean or the median value on the block. When there is an object, we need to know if it is darker or lighter than the background. In our cases, the bias is multiplicative, and we take a mean operator or a median one (in stead of a max or a min for an additive bias).

We obtain an image of size $N \times N$ used in an interpolation to reconstruct the total bias on each fiber. We used a bilinear interpolation. The *libirimage* library could be completed to implement a cubic spline interpolation.

We could note that the bias can be estimated more precisely using other means more subtle. Our method stay fast and seems sufficient.

Once the bias is estimated, we need to divide the image by its bias. We need to check that the bias is above one to avoid multiplication in fact ! In the calibration image we [...]

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will restore a false dynamic of 1000 levels to avoid any issues with the second division for the calibration of the fiber injection rate. (cf. pat 1.4.1).

1.4 Calibration and reconstruction bloc

1.4.1 Injection rate calibration

The flow estimation block on the left is used to have an image representing the calibration rate per fiber. We can use it to restore a 100% injection on all fibers. To do this, we divide the object image by the injection rate after some filtering on the calibration values (greater than 1, with enough dynamics).

1.4.2 Mosaic reconstruction

The mosaic reconstruction is obtained by spreading on all the fiber surface the same flow value estimated after the calibration image.

1.4.3 RBF reconstruction

To do ...

1.5 Results and discussion

The results are visible on the image database...

1.5.1

Figures 1.5 and 1.6 show all steps of the image processing. In this case, the background of the calibration image is estimated, that is why it does not appear on the images. The object background is given. The calibration images, and the object are unbiased. This results was obtained using default parameters. On this example : we note that the bias images have not the usual circular symmetry. The reason is that the main bias comes from the detector saturation resulting in a new additive background, not present in our first estimations, but the nature of which was attenuated by the processing.

1.5.2 Discussion

Numerous ideas of improvements have been proposed in this chapter. Some of them are more important than others. Here is the proposition of ideas by decreasing order of importance :

Reconstruction using Radial Basis Functions : this should results in smoother images, and allow for further processing

Remove the image zoom : the zoom is useful to lower the effect of the numerical opening before the watershed. We need to find the equivalent operator on an unzoomed image. An hypothesis to dig is that it would be equivalent with an opening using a cross-shape structuring element (?).

Fiber detection auto-control : to know if the detection is correct, or to check if the calibration need to be re-estimated on an acquisition.

Background auto-control : to know if the proposed background is correct, or if it needs to be re-estimated. In this case, the background is not given, do we need one ?

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(a) Initial image of the fibers (b) connected components of the fibers

(c) flow estimation on (b) - (d) subtraction estimated on (c)

(e) un-biased image from (d) - (f) background of the object (see fig. 1.6)

Fig 1.5 - The different steps of the imagecell process.

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(e) unbiased image from (d) - (f) final result

Fig 1.6 Different steps of the imagecell process.

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Continuous improvement of the fiber detection : we could use only the object images to improve or correct the fiber detection

Bias estimation using splines : the bias is estimated by bi-linear re-sampling that produces a bloc effect that could be embarrassing

We could also anticipate several research axis for future features : this will be addressed in an other report.